

# Colposcope E For the gynecologist's office.



# Colposcope E - Low price, high quality.

The compact Colposcope E from Carl Zeiss proves that a flexible, very convenient colposcope can also provide high optical quality.

And it also proves that "Economy" does not automatically mean less convenience. But it does mean an economical price.

The following models are available, all tailored precisely to your requirements in the gynecologist's office:

- On the suspension mount/on the floor stand
- With a straight tube/with an inclined tube
- With a fixed magnification/with a 3-step magnification changer
- Video compatible:

The MediLive® Primo 1CCD video camera (option) provides impressive information on the patient and superbly documents the results of the examination with the Colposcope E.

## Technical Data

**Illumination:** 6V, 15W halogen illumination

**Working distance:** 300 mm

**Magnification:** 5x or with 3-step magnification changer: 3x, 5x, 8x

**Green filter:** swing-in

**Video camera:** 1CCD, 1/2", S-Video/VBS, PAL/NTSC

**Suspension system:** suspension mount (connecting piece required for examination chair\*), mobile floor stand

\*Mounts for attaching connecting pieces to examination chairs are available from the chair manufacturer.

## Compliance:

- DIN EN ISO 9001
- DIN EN ISO 13485
- DIN EN 60601-1



*Application photos:  
Prof. Hillemanns, University Freiburg,  
Gynecological Hospital, Germany*

## Carl Zeiss Surgical GmbH

A Carl Zeiss Meditec Company  
73446 Oberkochen  
Fax: +49 (0) 7364/20-4823  
Email: [surgical@meditec.zeiss.com](mailto:surgical@meditec.zeiss.com)  
[www.meditec.zeiss.com/colposcopy](http://www.meditec.zeiss.com/colposcopy)  
[www.meditec.zeiss.com/contacts](http://www.meditec.zeiss.com/contacts)

MediLive is a registered trademark of Carl Zeiss.